

IV Administration Set

With Universal Spike and Rotating Luer

Our IV Administration Set with universal spike, male swivel luer lock, roller clamp, and drip chamber with 15u filter. Each set is individually blister-packaged in quantities of 50 per box and features patient information directly on the exterior of the package.



ACC-F001

- 20 drops/ml
- One (1) Y-injection site 6" from distal end
- Male swivel luer lock connector
- Universal vented or non-vented spike
- Roller clamp
- 15 micron filter located in drip chamber
- 92" length
- DEHP-free
- Latex-free

Manufactured by **TrueCare** USA 
BIOMEDIX



Access Medical USA • 10440 Belga Dr. San Antonio, Texas • 78240
www.accessmedicalusa.com

800-599-0409