

# IV Administration Set

With GVS Speedflow® .22u In Line Filter

Our IV Administration Set with a 0.2 Micron Filter consists of a universal spike, male swivel luer lock, roller clamp, drip chamber with a 15 micron filter, and a .02 micron in line filter. Each set is individually blister-packaged in quantities of 50 per box and features patient information directly on the exterior of the package.



## ACC-F022

- 20 drops/ml
- One (1) Y-injection site 6" from distal end
- Male swivel luer lock connector
- Universal vented or non-vented spike
- Roller clamp
- 15 micron filter located in drip chamber
- 0.2 micron GVS Speedflow® in line Filter
- 92" length
- DEHP-free
- Latex-free

Manufactured by **TrueCare** USA  
BIOMEDIX



Access Medical USA • 10440 Belga Dr. San Antonio, Texas • 78240  
www.accessmedicalusa.com

**800-599-0409**