


IV Administration Set

With GVS EasyDrop® Flow Regulator

Our IV Administration Set with Flow Regulator consists of a universal spike, male swivel luer lock, roller clamp, flow regulator, and drip chamber. Each set is individually blister-packaged in quantities of 40 per box and features patient information directly on the exterior of the package.



ACC-F033

- 20 drops/ml
- One (1) Y-injection site 6" from distal end
- Male swivel luer lock connector
- Universal vented or non-vented spike
- Roller clamp
- 15 micron filter located in drip chamber
- Smooth-moving flow regulator dial range of 5ml-250ml per hour
- 92" length
- DEHP-free
- Latex-free

Manufactured by
TrueCare USA
BIOMEDIX



Access Medical USA • 10440 Belga Dr. San Antonio, Texas • 78240
www.accessmedicalusa.com

800-599-0409